

Portable Calibrator Versa

Easy and efficient calibration of ergometers



Highlights

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

The Calibrator Ultra is an easy to set up and very accurate portable calibration station. The calibration station is controlled by a PC or laptop. With this station it is possible to control or calibrate all existing ergometers accurately. The Lode Calibrator is mounted on the axis of the ergometer. This system is approved by an independent institute, and is used by calibration institutes all over the world since 1993. Calibrations can be done up to 750 watt. Due to the possibilities you create when using a normal PC to calibrate the calibrator you can print a calibration report for your customer and/or create your own database of calibration data. The Calibrator Ultra is *the* essential tool for professionally and efficiently servicing ergometers.

The Calibrator Versa is the next development in professional portable calibration equipment, coming from a pedigree of 2 decades of calibration devices supporting accurate ergometry.

Please note, specifications at 110 V are different than at 230 V.



Portable Calibrator Versa

Easy and efficient calibration of ergometers



Features



Extremely Accurate Calibration of Ergometers

Highly accurate calibration of ergometers is possible with this product



USB connectivity

USB to connect to PC or ECG or ergospirometry products facilitates easy connectivity.

Professional
tool

Professional Service Tool

Lode provides tools for service engineers of distributors that are completely dedicated to its products and purpose. Some can even be used for non-Lode calibration activities

Universal

Universal Calibration

This tool can be used for Calibration of Lode products as well as products of other manufacturers

Multiple
connectors

Multiple ergometer axis connectors available

Various accessories are available to connect the calibrator to most available ergometers

Self control

With the included balance arm and weight basket it is possible to control the unit itself. The Calibrator 2000 is standard supplied with a toolbox with adaptors for calibrating all current Lode ergometers (excl. MRI ergometers).










The Calibrator is supplied with the latest edition of Windows-based software.

Portable Calibrator Versa

Easy and efficient calibration of ergometers



Portable Calibrator Versa can a.o be extended with the following options:

<p>Universal adapter for pinless crank</p> <p>Calibrating various ergometers</p>  <p>Partnumber: 919028</p>	<p>Universal adapter for 1/2" thread crank</p> <p>Calibrating various ergometers</p>  <p>Partnumber: 919805</p>	<p>Adapter for crank with pin</p> <p>Calibrating various ergometers</p>  <p>Partnumber: 919051</p>	<p>Adapter for Dynavit crank</p> <p>Calibrating various ergometers</p>  <p>Partnumber: 919803</p>	<p>USB to Serial converter</p> <p>Easy connection</p>  <p>Partnumber: 226012</p>
<p>Breaking pins, set of 15</p> <p>Additional spare breaking pins</p>  <p>Partnumber: 919121</p>	<p>Isis adapter m8x1</p> <p>Calibrating various ergometers</p>  <p>Partnumber: 919095</p>	<p>RS232 cable</p> <p>Easy connection</p>  <p>Partnumber: 930911</p>	<p>Transportcase Calibrator Ultra / Versa</p> <p>Safe transport</p>  <p>Partnumber: 919810</p>	

Portable Calibrator Versa

Easy and efficient calibration of ergometers



Specifications

Workload

Minimum rotational speed	25 rpm
Maximum rotational speed	180 rpm
Rotation direction clockwise	✓
Rotation direction anti-clockwise	✓
Maximum motor torque	55 Nm
Maximum sensor torque	100 Nm
Warming up time	5 min
Approved power range (maximum)	750 W
Approved power range (minimum)	0 W
Power range at maximum rpm (maximum)	750 W
Power range at maximum rpm (minimum)	0 W
Crank adapter included	✓

Accuracy

Speed accuracy	0.8 %
Power accuracy error > 100 W	1 %
Power accuracy error up to 100W	1 W
Torque resolution	0.05 Nm
Time resolving	10 ms
Accuracy rpm adjustment	1.0 rpm
Rpm resolution	1 rpm

User Interface

Software available in English	✓
Software available in German	✓
Software available in Dutch	✓
Software available in Japanese	✓

Connectivity

RS232 in connector	✓
--------------------	---

Dimensions

Cable length	3000 mm	118.1 inch
Product length (cm)	52 cm	20.5 inch
Product width (cm)	32 cm	12.6 inch
Product height	52 cm	20.5 inch
Product weight	29 kg	63.9 lbs

Power requirements

230 VAC	10 A
110 VAC	15 A
Phases	1
Frequency	50/60 Hz

Order info

Partnumber: 919940

*Specifications are subject to change without notice.

Standards & Safety

Integrated emergency stop ✓

Minimum System Requirements

Intel Core based processor	2.4 GHz
Monitor	1920 x 1080 pixels
Internal RAM	4000 MB
Minimum Free Hard Disk space	20 GB
Number of free USB ports	2
Mouse	✓
Keyboard	✓

Compatibility

Suitable for all brands	✓
Suitable for high power ergometers	✓
Microsoft Windows 10	✓

Included parts

PC software included ✓



FOR LIFE | UNDERSTANDING
MOVEMENT & PERFORMANCE

ISO 9001:2016 and ISO 13485:2016

Lode B.V.
Zernikepark 16
9747 AN Groningen
The Netherlands
Tel: +31 50 5712811
Fax: +31 50 5716746
E-mail: ask@lode.nl
Internet: www.lode.nl