

Heart rate for treadmills

Heart rate in beats per minute

Original
Quality
Option



Heart rate in beats per minute

This option gives the heart rate in beats per minute on the display of the control unit and offers the possibility of heart rate controlled walking on all treadmill models. The test person wears a belt with a sensor. Note: this is only possible in combination with the control unit option installed. The heart rate receiver is mounted on the handlebar

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



Original Accessory

Accessory designed and manufactured according to the Lode company quality standards

Specifications

Accuracy

Maximum heartrate	240 bpm
Minimum heartrate	15 bpm

Minimum System Requirements

Maximum distance from receiver	80 cm
--------------------------------	-------

Order info

Partnumber: 945801

**Specifications are subject to change without notice.*