


Adjustable sports cranks incl. pediatric range

Original
Quality
Option

Optimal force application



Optimal force application

Special step less adjustable cranks to be built-in in the Lode factory. Range 95 - 185 mm; minimum leg length: 650 mm, commercial available click pedals can be mounted. This option will be supplied instead of the standard cranks.

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



Original Accessory

Accessory designed and manufactured according to the Lode company quality standards

Specifications

Comfort		
Minimum leg length user	635 mm	25 inch
Crank adjustment length	90 mm	3.5 inch
Maximum crank length	185 mm	7.3 inch
Minimum crank length	95 mm	3.7 inch

Order info

Partnumber: 925808

**Specifications are subject to change without notice.*



ISO 9001:2016 and ISO 13485:2016

Lode B.V.
Zernikepark 16
9747 AN Groningen
The Netherlands
Tel: +31 50 5712811
Fax: +31 50 5716746
E-mail: ask@lode.nl
Internet: www.lode.nl