


# Pedal shoes (pair)


Extra stability during cycling

Original  
Quality  
Option



## Extra stability during cycling

A pair of pedal shoes for foot fixation to enable reclined pedaling. These pedal shoes will replace the standard pedals of the ergometer. If you would like to order them as additional to the standard pedal, please use part number P917803 on your order.



## Multifunctional

The ergometer can be used in various ergometry settings, enabling a multifunctional deployment.

## Stability

Since athletes are getting more and more powerful and testing more extreme, the ergometer is designed for high workloads up to 2500 watt. Even the strongest and most powerful athletes will experience a stable and comfortable basis.

**7**  
**watt**

### Extreme low start up load

The extreme low start-up load of 7 watts and the adjustability in small steps of 1 watt make this ergometer perfectly suitable for many different applications. The standard control unit shows multiple ergometry parameters and you can determine your specific default setting and start-up menu.



# Pedal shoes (pair)

Extra stability during cycling

Original  
Quality  
Option

## Specifications

### Comfort

Shoe size pedal shoe EU	32 - 41
Shoe size pedal shoe US male	1 - 8
Shoe size pedal shoe US female	1 - 9
Shoe size pedal shoe UK	12.5 (ch) - 7.5 (ad)

### Dimensions

Product length (mm)	240 mm
Product width (mm)	90 mm

### Order info

Partnumber: 917803

*\*Specifications are subject to change without notice.*