


Pedal shoes pediatric (pair)


Pedal shoes for children

Original
Quality
Option



Pedal shoes for children

Pedal shoes allowing children to do reclining ergometry



Multifunctional

The ergometer can be used in various ergometry settings, enabling a multifunctional deployment.

Stability

Since athletes are getting more and more powerful and testing more extreme, the ergometer is designed for high workloads up to 2500 watt. Even the strongest and most powerful athletes will experience a stable and comfortable basis.

**7
watt**

Extreme low start up load

The extreme low start-up load of 7 watts and the adjustability in small steps of 1 watt make this ergometer perfectly suitable for many different applications. The standard control unit shows multiple ergometry parameters and you can determine your specific default setting and start-up menu.



Pedal shoes pediatric (pair)

Pedal shoes for children

Original
Quality
Option

Specifications

Comfort

Shoe size pedal shoe EU	24 - 31
Shoe size pedal shoe US male	8 - 13
Shoe size pedal shoe UK	6.5 - 12

Dimensions

Product length (mm)	170 mm
Product width (mm)	60 mm

Order info

Partnumber: 917833

**Specifications are subject to change without notice.*



FOR LIFE | UNDERSTANDING
MOVEMENT & PERFORMANCE

ISO 9001:2016 and ISO 13485:2016

Lode B.V.
Zernikepark 16
9747 AN Groningen
The Netherlands
Tel: +31 50 5712811
Fax: +31 50 5716746
E-mail: ask@lode.nl
Internet: www.lode.nl