

# Remote Control for treadmill

Ease of use

Original  
Quality  
Option



## Ease of use

Remote Control for Lode treadmills. It can be used for controlling the speed and elevation as well as for starting and stopping. Note: can only be ordered ex factory in addition to 945814 or 945815.

## Easy to operate

For Lode products this means:

- easy to connect
- easy to move around
- easy user interface

## High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

## Multifunctional

The treadmill can be used in various ergometry fields, like cardiology, pulmonary function, rehabilitation, sports medicine and physiotherapy, as a stand alone unit or controlled by other equipment.



## Original Accessory

Accessory designed and manufactured according to the Lode company quality standards



# Remote Control for treadmill

Ease of use

Original  
Quality  
Option

## Specifications

### Dimensions

Cable length	2990 mm	117.7 inch
Product height (mm)	25 mm	
Product length (mm)	100 mm	
Product width (mm)	47 mm	
Product weight	0.210 kg	0.5 lbs

### Order info

Partnumber: 945831

*\*Specifications are subject to change without notice.*