

# Shortened saddle shaft

Increase flexibility for smaller people

Original  
Quality  
Option



## Increase flexibility for smaller people

With the shortened saddle shaft, the saddle of the ergometer can be adjusted to the lowest position possible. The saddle is adjustable in a range of 452 - 775 mm. The minimum leg length (with standard pedals) is 622 mm and with adjustable cranks 532 mm. The minimum saddle height is 452 mm. NB: this option cannot be combined with the electric adjustable saddle height (#960810)



## Multifunctional

The ergometer can be used in various ergometry settings, enabling a multifunctional deployment.

## High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



## Original Accessory

Accessory designed and manufactured according to the Lode company quality standards



# Shortened saddle shaft

Increase flexibility for smaller people

Original  
Quality  
Option

## Specifications

### Comfort

Minimum leg length user (incl. adjustable pedals)

532 mm

20.9 inch

### Order info

Partnumber: 960806

*\*Specifications are subject to change without notice.*



**FOR LIFE** | UNDERSTANDING  
MOVEMENT & PERFORMANCE

ISO 9001:2016 and ISO 13485:2016

Lode B.V.  
Zernikepark 16  
9747 AN Groningen  
The Netherlands  
Tel: +31 50 5712811  
Fax: +31 50 5716746  
E-mail: ask@lode.nl  
Internet: www.lode.nl