

Shortened saddle shaft

Increase flexibility for smaller people

Original
Quality
Option



Increase flexibility for smaller people

With the shortened saddle shaft, the saddle of the ergometer can be adjusted to the lowest position possible. The saddle is adjustable in a range of 452 - 775 mm. The minimum leg length (with standard pedals) is 622 mm and with adjustable cranks 532 mm. The minimum saddle height is 452 mm. NB: this option cannot be combined with the electric adjustable saddle height (#960810)



Multifunctional

The ergometer can be used in various ergometry settings, enabling a multifunctional deployment.

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



Original Accessory

Accessory designed and manufactured according to the Lode company quality standards



Shortened saddle shaft

Increase flexibility for smaller people

Original
Quality
Option

Specifications

Comfort

Minimum leg length user (incl. adjustable pedals)

532 mm

20.9 inch

Order info

Partnumber: 960806

**Specifications are subject to change without notice.*



FOR LIFE | UNDERSTANDING
MOVEMENT & PERFORMANCE

ISO 9001:2016 and ISO 13485:2016

Lode B.V.
Zernikepark 16
9747 AN Groningen
The Netherlands
Tel: +31 50 5712811
Fax: +31 50 5716746
E-mail: ask@lode.nl
Internet: www.lode.nl