

Adjustable Acceleration & Deceleration for Katana & Ultra

Original
Quality
Option

Flexibility in exercising



Flexibility in exercising

This option makes it possible to adjust the speed of the acceleration (and deceleration) of the Katana and Ultra treadmill. The standard acceleration is set to 6 seconds per 10 km. This option makes it possible to set the value between 3 and 24 seconds per 10 km.

Smooth acceleration

All treadmills have a smooth acceleration from 0 km/h to the start speed.

Multifunctional

The treadmill can be used in various ergometry fields, like cardiology, pulmonary function, rehabilitation, sports medicine and physiotherapy, as a stand alone unit or controlled by other equipment.

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



Low noise

Material choices, refined components and accurate manufacturing techniques lead to low noise.

0.1
km/h

Small adjustment steps

The speed of the treadmill can be adjusted in the smallest steps of only 0.1 km/h!



Adjustable Acceleration & Deceleration for Katana & Ultra

Original
Quality
Option

Flexibility in exercising

Speed adjustment steps

0.1 km/h

0.1 mph

Order info

Partnumber: 945847

**Specifications are subject to change without notice.*



FOR LIFE | UNDERSTANDING
MOVEMENT & PERFORMANCE

ISO 9001:2016 and ISO 13485:2016

Lode B.V.
Zernikepark 16
9747 AN Groningen
The Netherlands
Tel: +31 50 5712811
Fax: +31 50 5716746
E-mail: ask@lode.nl
Internet: www.lode.nl