

# Bluetooth Smart heart rate

Heartrate available within an extreme wide range

Original  
Quality  
Option



## Heartrate available within an extreme wide range

This option gives the heart rate in beats per minute on the display of the control unit and offers the possibility of heart rate controlled testing or exercising on all treadmill models and ergometers. Thanks to the bluetooth, the heart rate can be received over a long distance, up till 10 meters, without interference. The test person wears a belt with a sensor. It can only be ordered in combination with the (Programmable) control units with touchscreen (945814/5 and 945834/5). This option includes 1 belt. Please note that this option is NOT suitable for triggering of Blood Pressure measurement. If you want to measure the HR in combination with the BP, you can use the HR signal of the custom pack (#945888). The Bluetooth HR option is not possible then.

### High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

### Easy to operate

For Lode products this means:

- easy to connect
- easy to move around
- easy user interface

### Freedom of Movement

The test subject can move freely without any obstruction of the system



### Original Accessory

Accessory designed and manufactured according to the Lode company quality standards

### Bluetooth connectivity



Bluetooth connectivity allows for secure short and medium range connectivity of peripherals



# Bluetooth Smart heart rate

Heartrate available within an extreme wide range

Original  
Quality  
Option

## Specifications

### Accuracy

Maximum heartrate 240 bpm

Minimum heartrate 15 bpm

### Minimum System Requirements

Maximum distance from receiver 1000 cm

### Order info

Partnumber: 945833

*\*Specifications are subject to change without notice.*



**FOR LIFE** | UNDERSTANDING  
MOVEMENT & PERFORMANCE

ISO 9001:2016 and ISO 13485:2016

Lode B.V.  
Zernikepark 16  
9747 AN Groningen  
The Netherlands  
Tel: +31 50 5712811  
Fax: +31 50 5716746  
E-mail: ask@lode.nl  
Internet: www.lode.nl