

Adjustable cranks

Optimal force application

Original
Quality
Option



Optimal force application

Adjustability of cranks to allow for optimal force application at different limb length. The cranks are adjustable in a range of 80 - 180 mm, in steps of 5 mm. This option will be supplied instead of the standard cranks.



Multifunctional

The ergometer can be used in various ergometry settings, enabling a multifunctional deployment.

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2015, and ISO 13485:2016 certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.



Original Accessory

Accessory designed and manufactured according to the Lode company quality standards



Adjustable cranks

Optimal force application

Original
Quality
Option

Specifications

Comfort

Crank adjustment steps	5 mm	0.2 inch
Crank adjustment length	100 mm	3.9 inch
Maximum crank length	180 mm	7.1 inch
Minimum crank length	80 mm	3.1 inch

Order info

Partnumber: 960804

**Specifications are subject to change without notice.*