



Ref. Certif. No.

HU-002756

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB SCHEME

CB TEST CERTIFICATE

Product

Cycling ergometer with Lode Ergometry Manager (LEM) software

Name and address of the applicant

LODE B.V.  
Zernikepark 16, 9747 AN Groningen,  
The Netherlands

Name and address of the manufacturer

LODE B.V.  
Zernikepark 16, 9747 AN Groningen,  
The Netherlands

Name and address of the factory

LODE B.V.  
Zernikepark 16, 9747 AN Groningen,  
The Netherlands

Note: When more than one factory, please report on page 2

Ratings and principal characteristics

100-240V~, 50/60Hz; 160W;  
Class I, Type B applied part.  
(see Test Report for further ratings)

Trademark (if any)



Customer's Testing Facility (CTF) Stage used

---

Model / Type Ref.

Excalibur series cycling ergometers:  
- Excalibur CPET (965900), Brachumera CPET (965901),  
- Excalibur Sport (965910), Brachumera Sport (965911),  
- Excalibur Sport PFM (965920),  
Brachumera Sport PFM (965921)  
Software:- Lode Ergometry Manager (LEM) software (955900)  
(see Test Report for further details)

Additional information (if necessary may also be reported on page 2)

---

A sample of the product was tested and found to be in conformity with

IEC 60601-1:2005 + COR1:2006 + COR2:2007 + AMD1:2012  
IEC 60601-1-6:2010 + AMD1:2013  
National Differences: CA, US

As shown in the Test Report Ref. No. which forms part of this Certificate

28250166 001

This CB Test Certificate is issued by the National Certification Body



TÜV Rheinland InterCert Kft., Division MEEI  
H-1132 Budapest, Váci út 48/A-B  
www.tuv.com

Date: 2019-10-25

Signature:

*Gabor Kassai*  
Gabor Kassai