

IEC**IECEE
CB
SCHEME**

Ref. Certif. No.

HU-002105

**IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT
(IECEE) CB SCHEME****CB TEST CERTIFICATE**

Product

Ergometer family

Name and address of the applicant

LODE B.V.
Zernikepark 16, 9747 AN Groningen,
The Netherlands

Name and address of the manufacturer

LODE B.V.
Zernikepark 16, 9747 AN Groningen,
The Netherlands

Name and address of the factory

LODE B.V.
Zernikepark 16, 9747 AN Groningen,
The Netherlands

Note: When more than one factory, please report on page 2

Ratings and principal characteristics

100-240V~; 50/60Hz; 160W; Class I;
Type B applied part

Trademark (if any))



Customer's Testing Facility (CTF) Stage used

Model / Type Ref.

Corival Recumbent series:
969900 - CPET; 969903 - Pediatric;
969901 - Rehab; 969912 - OEM;
969902 - Flow; 969904, 969906 - Connect;

Additional information (if necessary may also be reported on page 2)

A sample of the product was tested and found to be in conformity with

IEC 60601-1:2005 + COR1:2006 + COR2:2007 + AMD1:2012
IEC 60601-1-6:2010 + AMD1:2013
National Differences: CA, US

As shown in the Test Report Ref. No. which forms part of this Certificate

28241612 001

This CB Test Certificate is issued by the National Certification Body

TÜV Rheinland InterCert Kft., Division MEE
H-1132 Budapest, Váci út 48/A-B
www.tuv.hu

Date: 2017-12-01

Signature:

Gabor Kassai